



EXPORT PRODUCT CATALOGUE 2024



MDT (est. 1999) is a reliable manufacturer and supplier of ophthalmic medical devices, accessories and consumables (including BioGlo™ Fluorescein Sodium Strips, TearFlo™ Schirmer Test Strips, and CRP - Chin Rest Papers).

Thanks to our experience, a long-lasting international presence, sustainable development, focus on craftsmanship and excellent business partners we are able to meet the growing requirements by offering a wide range of high-quality, functional products for medical practitioners and OEM market.

Among the others we manufacture ophthalmic surgical tables, operator's chairs, refraction units (workstations), electric tables, stools, as well as devices used for diagnostics and treatment of strabismus and amblyopia, which are proudly 'made in EU' and CE marked, as well as FDA approved.

We are proud of our Authorised Service Partners and service engineers who provide the customers with an outstanding level of technical support, both at the stage of equipment installation and in case of any issues which could have occurred during a product life-time.

Should you need any further information, please do not hesitate to contact your local business partner, as they look forward to hearing from you!

MDT Team





Simplicity and space-saving design combined with functionality.

TECH SPEC


- Equipped with electrically elevated SLIM chair (available in black upholstery only)
- Convenient chair height adjustment switch on the table top
- V-shape, fixed height, non-extendable table top for 2 instruments, with optional universal chin rest adapter
- ALU post suitable for installation of both, chart projector and PA-1000 phoropter arm (option)
- Table top electromagnet lock switch hidden under the auxiliary table top
- Efficient, foldable and revolving LED lamp with a manual switch
- Right-handed setup as a standard, left-handed available upon order or during on-site installation
- Side casings and table top painted in RAL7047
- Front casing and auxiliary table top painted in RAL7024
- Steel base powder coated in RAL7024 EU Declaration of Conformity provided



PA-1000




UNIQ with SLIM chair (right-handed setup)



Left-handed setup*



Chair height adjustment switch



Pivotal arm with trial lenses tray and cassette**/**



LED lamp switch



Table top electromagnet lock switch



C-13 sockets



CBS-1000

*- optional feature
**- trial lens set – optional feature

OPTIONAL EXTRAS

In order to provide the most competitive prices the number of optional features has been limited:

- Footrest for SLIM chair - n/a
- COMFORT chair – n/a
- PA-1000 counterbalanced phoropter arm equipped with a gas spring
- PA-3000 - n/a
- Elevated table top – n/a
- Customised main body colours - n/a
- Customised upholstery colours - n/a
- Pivotal arm with trial lenses tray and cassette
- Chart projector shelf
- Left-handed setup (factory preset or upon on-site installation)
- CBS-1000 Control Box Stand
- Ophthalmoscope/retinoscope rack
- Universal chin rest adapter

Dimensions	SLIM CHAIR
Width	1 145mm
Depth	1 100 - 1 610mm (±10)
Height	1 630mm
Chair height working range	430 - 630mm (±5)
Table top dimensions	W 944/900/735mm x D 445/410/445mm
Table top height (fixed)	861mm
Table top working range (height adjustable)	n/a

Technical specification

Supply voltage	100 - 240V AC, 50/60Hz
Power	280VA
LED lamp power	9W
Power supply cord	3m with UniSchuko plug
Max. table top load	50kg
Max. chair load	150kg
Max. power of devices connected to 2x C-13 sockets	450VA
Packaging	1 pallet, 120x80x160cm, ca. 150kg (gross)

UNIQ WITH SLIM CHAIR

