



## The OCULUS LidStick® and Dispenser

The LidStick® was designed to replace the abrasive and uncomfortable cotton and wood applicators. Utilizing medical grade silicone and plastic to become the ultimate multi-purpose tool for your eye care practice.

The LidStick® was developed by OCULUS to manipulate the eyelids for dry eye measurements, slit lamp examinations, cornea-scleral surface detection and more.

One common tool for 5 key functions:

- meibography
- cornea-scleral scans
- sensitive squeezing of the meibomian glands
- careful debriding of the lid margins
- estimate the sizes of a chalazion etc.



For more information please watch this video:

<https://www.youtube.com/watch?v=tT6fx2RSQis&pp=ygUPb2N1bHVzIGxpZHN0aWNr>



The Dispenser is a practical storage container for hygienic individual removal of the LidStick®. It combines the storage of the LidStick® with the function of removing them hygienically. The specially designed aluminum storage container holds the capacity of one roll of LidStick®.

The guide rail enables the individual LidStick® to be torn off in a targeted manner and prevents uncontrolled winding or unwinding of the roll. With a viewing window, the quantity of individual LidStick® can be easily estimated from the side caps.

Whether lying, standing or hanging - the dispenser can be easily attached with the enclosed PosterStrips.

