

49700



48110



## Retinoscopy Racks/ -Mirror

**49700** Retinoscopy mirror, Ø 45 mm

**48110** Retinoscopy racks with slider (by Trousseau), in frame of light metal and strengthened synthetic plastic material, set consists of 2 racks, plus and minus from 0.5D - 10D in half diopter steps, from 10 - 20D in 1 D steps, in case

36620



## Sphere-balance with handle, coated

**36610** sph +/- 0.12 dpt

**36620** sph +/- 0.25 dpt

**36640** sph +/- 0.50 dpt

36600



## Balance Rack

**36600** Sphere-balance rack, sph +/- 0.25 and +/- 0.50 dpt

36730



## Vario-Cross-Cylinder

**36730** Vario-cross-cylinder (by Reiner), with effect changes +/- 0.25 dpt and +/- 0.50 dpt, special coated

36720



## Cross Cylinder, with handle, AR-coated

**36710** +/- 0.12 dpt

**36720** +/- 0.25 dpt

**36740** +/- 0.50 dpt

**36780** +/- 1.00 dpt